



The RNA-seq Protocol

Erin Osborne Nishimura
DSCI 512:
RNA-seq data analysis
November 5, 2024

The evolving face of High Throughput Sequencing

- **Short reads** – Illumina (aka **nextgen** sequencing)
- **Third gen sequencing** -- Long reads -- Pac Bio & minION (Oxford Nanopore Technologies)

illumina®




PACBIO®

Oxford
NANOPORE
Technologies®

Genomics

Three flavors of high-throughput sequencing





illumina®

pb PACBIO®


Oxford
NANOPORE
Technologies®

Pac Bio – circular long-read sequencing

- Immobilized DNA Polymerases at the bottom of each well amplify a single circularized DNA molecule.
- Colored light is emitted for each nucleotide
- Long reads are produced
- Errors are minimized by going around each circle multiple times



Explore Pac Bio Sequencing on your own



https://youtu.be/_ID8JyAbwEo?si=-4pI3dLq5AoMlmeZ


<https://www.pacb.com/blog/long-read-sequencing/>

Oxford Nanopore – MinION




Oxford Nanopore sequencing technology, a third-generation sequencing technology

- minION
- A single, long, unamplified DNA molecule is dissociated into a single strand and fed through a small protein pore in a membrane.
- As the DNA passes through the membrane, a small current is created that is distinct for each nucleotide.
- The ionic flow signature of each nucleotide can be read through an extremely sensitive detector




Nanopore sequencing technology



Genomics

Three flavors of high-throughput sequencing




























illumina®

pb PACBIO®


Oxford
NANOPORE
Technologies®

Illumina Sequencing



	Benchtop Sequencers		Pico- & Scale Sequencers		
					
	iSeq 100	MiniSeq	MiSeq Series +	NextSeq 550 Series +	NextSeq 1000 & 2000
Popular Applications & Methods	Key Application 	Key Application 	Key Application 	Key Application 	Key Application 
Large Whole-Genome Sequencing (human, plant, animal)					
Small Whole-Genome Sequencing (microbe, virus)					
Exome & Large Panel Sequencing (enrichment-based)					
Targeted Gene Sequencing (amplicon-based, gene panel)					
Single-Cell Profiling (scRNA-Seq, scDNA-Seq, oligo tagging assays)					
Run Time	9.5–19 hrs	4–24 hours	4–55 hours	12–30 hours	11–48 hours
Maximum Output	1.2 Gb	7.5 Gb	15 Gb	120 Gb	360 Gb*
Maximum Reads Per Run	4 million	25 million	25 million †	400 million	1.2 billion*
Maximum Read Length	2 × 150 bp	2 × 150 bp	2 × 300 bp	2 × 150 bp	2 × 150 bp

Illumina Sequencing



Benchtop Sequencers

Production-Scale Sequencers



School of Medicine
UNIVERSITY OF COLORADO
ANSCHUTZ MEDICAL CAMPUS



NextSeq 550 Series



NextSeq 1000 & 2000



NovaSeq 6000

Popular Applications & Methods	Key Application	Key Application	Key Application
Large Whole-Genome Sequencing (human, plant, animal)			
Small Whole-Genome Sequencing (microbe, virus)			
Exome & Large Panel Sequencing (enrichment-based)			
Targeted Gene Sequencing (amplicon-based, gene panel)			
Single-Cell Profiling (scRNA-Seq, scDNA-Seq, oligo tagging assays)			
Transcriptome Sequencing (total RNA-Seq, mRNA-Seq, gene expression profiling)			
Run Time	12–30 hours	11–48 hours	~13 – 38 hours (dual SP flow cells) ~13–25 hours (dual S1 flow cells) ~16–36 hours (dual S2 flow cells) ~44 hours (dual S4 flow cells)
Maximum Output			
Maximum Reads Per Run	120 Gb	360 Gb*	6000 Gb
Maximum Read Length	400 million	1.2 billion*	20 billion
	2 × 150 bp	2 × 150 bp	2 × 250**

Illumina short read sequencing: The basic protocol

RNA-seq library prep



1. total RNA




1. mRNA



1. cDNA



1. Library prep



2. Cluster generation



3. Sequencing by Synthesis



ATG CGGTTA TAGTAT A
TGC CACCTA ATTTC C
CAG GTCCCC TTTGTT A
ATG CGGTTA TAGTAT A
TGC CACCTA ATTTC C

4. Data analysis

- Total RNA extracted
- mRNA (2% of total) is enriched
- RNA converted to cDNA using Reverse Transcriptase
- NextGen sequenced
- Sequences are aligned to the genome. Exons will be represented.
- RNA processing and library prep are typically done with kits (NEB Next RNA-seq kit)
- NEXTGEN sequencing is often done at a core facility

Basic steps of Illumina sequencing

1. Library preparation

- fragment the DNA
- ligate adapters onto fragments




2. Cluster generation

- adhere fragments to a flow cell
- amplify them



3. Sequencing by synthesis


- directly image a specialized polymerization reaction that allows visualization of each nucleotide incorporated
- Tens of millions of fragments are sequenced




4. Data analysis

- Use computer algorithms to align sequences to the genome

ATGCGGTTATAGTATA
TGCCACCTAATTCC
CAGGTCCCCTTTGTTA
ATGCGGTTATAGTATA
TGCCACCTAATTCC




Step 1: Library preparation



- Total RNA is extracted
- mRNA is enriched either by mRNA enrichment or rRNA subtraction
 - Oligo dT beads
 - Ribozero, Ribominus, Ribogone
- mRNA is fragmented through heat + ions
- cDNA is produced by random oligo or oligo dT priming

Step 1 – Adding adapters

Adapters are ligated to the cDNA pieces. These are like handles that allow downstream instrumentation to interact with each molecule.



Please watch this video on your own
for homework:

- <https://www.neb.com/en-us/tools-and-resources/video-library?device=modal&videoid=%7Bd824c8c5-7942-437c-9086-e93ff3c94a12%7D>
- A link is on your homework assignment page

Step 1 - Multiplexing:

Multiple samples can be sequenced on the same flow cell







- Each sample is ligated to a specific adapter pair with its own sequence
- Samples are then merged and run on the same flow cell
- Samples are split computationally using unique adapter sequences

Adapters: A closer look


Unique dual index




-  **Flow cell binding sequence:** Platform-specific sequences for library binding to instrument
-  **Sequencing primer sites:** Binding sites for general sequencing primers
-  **Sample indexes:** Short sequences specific to a given sample library
-  **Insert:** Target DNA or RNA fragment from a given sample library

The aim of the sample prep step is to obtain nucleic acid fragments with adapters attached on both ends

RNA-seq



- Total RNA extracted
- mRNA (2% of total) is enriched
- RNA converted to cDNA using Reverse Transcriptase
- NextGen sequenced
- Sequences are aligned to the genome. Exons will be represented.
- RNA processing and library prep are typically done with kits (NEB Next RNA-seq kit)
- NEXTGEN sequencing is often done at a core facility




Illumina Sequencing: Step 2


2. Cluster generation

- adhere fragments to a flow cell
- amplify them


2 fragment capture and amplification

single-stranded fragments
bound to inside
surface of flow cell






Illumina Sequencing: Step 2



Illumina Sequencing: Step 2



iSeq 100



MiniSeq



MiSeq Series +



NextSeq 550 Series +



NextSeq 1000 & 2000



NextSeq 550 Series +



NextSeq 1000 & 2000



NovaSeq 6000





Illumina Sequencing: Step 3

3. Sequencing by synthesis

- Sequence using reversible dye termination.
- directly image the sequencing reaction
- Tens of millions of fragments are sequenced


3 sequencing





Illumina Sequencing: Step 3

Reversible dye termination incorporates one fluorescently labeled nucleotide at a time






Illumina Sequencing: Step 3


Reversible dye terminators

A




A regular deoxynucleotide (dNTP)

C



A dye-labeled dideoxynucleotide (ddNTP) used in Sanger sequencing

D



removable dye

A reversible dye terminator used in Sequencing By Synthesis


**removable
blocking
group**



Illumina Sequencing: Step 3

A picture of the flow cell is taken after the incorporation of each nucleotide

Base calling from raw data



What do we get out?



Single-end sequencing

- 150 bp sequenced on one side of the insert
- Depending on the organism & project, 10,000,000 – 200,000,000 inserts sequenced
 - 100 MB – 1 GB text file




Paired-end sequencing

- 150 bp sequence on each side of the insert
- Depending on the organism & the project, 10,000,000 – 200,000,000 inserts sequenced
 - Two 100 MB – 1 GB text files


4. Data analysis

Use computer algorithms to align sequences to the genome



- Sequenced fragments (called reads) will align to the exons
- There will be gaps at introns
- The number of reads mapped to a gene is proportional to its expression level (provided that some normalization has been performed)


RNA-seq data is visualized on a browser



- Genome is below. Gene models are below (3 genes are shown).
- Samples are rows (two wild-types and two mutants are shown).
- Track heights represent the number of aligned reads for that region.
- Click around on a web-interface for the whole genome


Question

A




A regular deoxynucleotide (dNTP)

C



A dye-labeled dideoxynucleotide (ddNTP) used in Sanger sequencing

D



A reversible dye terminator used in Sequencing By Synthesis

- What would happen if the sequencing tech accidentally substituted dNTPs for reversible dye terminators?
- What would happen if Sanger sequencing dye-labeled ddNTPs were used?

Decisions, Decisions



- Amount of input RNA?
- Method of mRNA enrichment?
- Strand specific (directional)?
- Type of library prep kit?
- Type of sequencing
 - Single-end or paired-end?
 - Read length, aka **cycles** (2x150)?
 - How many reads (> 20 million)?
 - What platform (novaSEQ or nextSEQ)?

Videos

- NEB Next Ultra Directional Library Prep
 - <https://www.neb.com/tools-and-resources/video-library?device=modal&videoid=%7Bd824c8c5-7942-437c-9086-e93ff3c94a12%7D>
- Illumina Sequencing By Synthesis
 - <https://www.youtube.com/watch?v=fCd6B5HRaZ8>
- If you're interested - 10x Genomics single-cell RNA-seq:
 - <https://www.youtube.com/watch?v=6UVOdCc1Q7I>